

## AF3113T Dino-Lite Basic Wireless




Dino-Lite **AF3113T Grey / White**, dirancang untuk memenuhi sebagian besar aplikasi dalam perbesaran 10x - 50x, 230x dilengkapi dengan kontrol kecerahan LED. Desain modular-nirkabel yang baru memungkinkan koneksi USB atau penggantian secara cepat ke WiFi Handphone dengan adaptor **WF-20 (Optional)** menjadikannya mikroskop tanpa kabel dan portabel. AF3113T tidak mempunyai fitur polarisasi tetapi tetap didukung software profesional **DinoCapture** (free update) untuk menyimpan gambar, video, dan memberikan keterangan sertamelakukan pengukuran.

**Series AF3113T adalah pengganti series AM3113T (Discontinue).**

Pertanyaan **contact support (Klik)**

Optional : Basic Stand & WF20



|                                  |                 |
|----------------------------------|-----------------|
| <b>Spesial feature</b>           | No              |
| <b>ESD Safe</b>                  | No              |
| <b>Measurement functionality</b> | Yes             |
| <b>Housing</b>                   | Regular         |
| <b>Cap Changeable</b>            | No              |
| <b>Polarizer</b>                 | No              |
| <b>Number of LEDs</b>            | 8               |
| <b>Light / LED type</b>          | White           |
| <b>Working distance</b>          | Standart        |
| <b>Resolution</b>                | VGA (640 x 480) |
| <b>Connection</b>                | USB 2.0         |

## ADDITIONAL INFO

|                            |   |
|----------------------------|---|
| Model                      | AF3113T Dino-Lite Premier   |
| Interface                  | USB 2.0   |
| Product Resolution         | VGA (640x480)   |
| Magnification              | 20x~55x, 220x   |
| Frame Rate                 | 30fps in VGA, YUY2<br>60fps in VGA, MJPEG   |
| Saving Format              | Image<br>DinoCapture2.0: BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM<br>DinoXscope: PNG, JPEG |
|                            | Video<br>DinoCapture2.0: WMV, FLV, SWF<br>DinoXscope: MOV   |
| Lighting                   | 8 white LEDs  |
| Microtouch                 | Yes   |
| Operating System Supported | Windows XP, Vista, 7, 8, 10<br>Mac OS 10.4 or later   |
| Calibration Function       | Yes   |
| Measurement Function       | Yes   |
| Unit Dimension             | 11.5cm (H) x 3.2cm (D)  |
| Unit Weight                | 105g  |
| Color                      | Grey: AF3113T / White: AF3113T(W)   |

