

NAIL

AM4113-N5UT (NAIL SCOPE)



Amati dan sirkulasi darah dokumen kuku untuk memeriksa kondisi arteri darah yang sehat.





SPECIFICATIONS

ADDITIONAL INFO

SPECIFICATIONS

Model	AM4113-N5UT Dino-Lite Premier
Interface	USB 2.0
Product Resolution	1.3M pixels (SXGA)
Magnification Rate	500x
Sensor	Color CMOS
Frame Rate	Up to 30fps
Image:	DinoCapture2.0: BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM DinoXscope: PNG, JPEG

Save Formats

Movie:
DinoCapture2.0: WMV, FLV, SWF
DinoXscope: MOV

Microtouch	Touch sensitive trigger on the microscope for taking pictures
Lighting	8 white LED lights switched on/off by software
Operating System Supported	Windows 8, 7, Vista, XP MAC OS 10.4 or later
Unit Weight	100(g)
Unit Dimension	9.4cm (H) x 3.2cm (D)
Package Dimensions	16cm (L) x 16cm (W) x 6cm (H)

I

Information about working distance and field of view

