

DINO EYE USB

AM7023B Dino-Eye (USB) (Mounting 30mm dan 30.5mm)



Dino-Eye Premier AM7023B adalah Eye Piece dengan interface USB 2.0 dan resolusi lima megapixel 2592x1944. Lensa ini cocok digabungkan dengan mikroskop tradisional dengan **slot lensa Diameter 30mm dan 30.5mm**.

Dino-Eye dapat digunakan di tempat produksi atau laboratorium untuk pemeriksaan dan dokumentasi.

Pertanyaan **contact support (Klik)**

Cara memilih KLIK

Warranty : 1 Tahun, Made in : Taiwan



SPECIFICATIONS

Image Quality	5.0 Megapixel
Connection	USB 2.0

ADDITIONAL INFO

SPECIFICATIONS

Model AM7023B Dino-Eye (USB)

Interface USB 2.0
Product Resolution 5M pixels (MJPEG codec may be required to run at 2592 x 1944 resolutions)
Sensor Color CMOS
Frame Rate Up to 30fps
Image:
DinoCapture2.0: BMP ,GIF ,PNG ,MNG ,TIF ,TGA ,PCX ,WBMP ,JP2 ,JPC ,JPG ,PGX ,RAS ,PNM

Save Formats

Movie:
DinoCapture2.0: WMV, FLV ,SWF

Measurement Function Yes

Calibration Function Yes

Operating System Supported Windows 7 ,Vista ,XP

Unit Weight 95(g) +/- 10(g)

Unit Dimension 5cm (H) x 3.4cm (D)

Fits 30mm/30.5mm diameter eyepiece slots

Package Dimensions 16cm (L) x 16cm (W) x 6cm (H)

