

AM7025X Dino-Eye Edge (USB)



Dino-Eyepieces Edge AM7025X mengadopsi kelas optik dengan sensor resolusi tinggi, dengan interface USB 2.0 dan resolusi 5 megapixel menawarkan performa gambar yang luar biasa serta menyediakan sudut pandang bidang relatif besar dari pendahulunya.

Dengan reproduktifitas warna yang halus, Dino-Eyepieces Edge AM7025X ini mampu menangkap detail terbaik dari gambar mikroskop, serta dapat digunakan dalam lingkup tempat produksi, atau laboratorium untuk pemeriksaan, dokumentasi, serta diterapkan pada aplikasi profesional.

Catatan: AM7025X kompatibel dengan DinoXcope versi 1,16, atau lebih baru

Dino-Eye Pertanyaan contact support (Klik)

Cara memilih KLIK

Warranty : 1 Tahun, Made in : Taiwan





SPECIFICATIONS

Image Quality	5.0 Megapixel
Diameter Adapter	23 mm,30mm,30.5mm
Connection	USB 2.0

ADDITIONAL INFO

SPECIFICATIONS

Model	AM7025X Dino-Eye Edge
Interface	USB 2.0
Product Resolution	5M pixels (2592x1944)
Frame Rate	10fps in 5MP/3MP/2MP, MJPEG 25fps in 1.3MP, MJPEG 30fps in VGA, MJPEG
Operating System Supported	Windows 10, 8, 7, Vista, XP Mac OS 10.9 or later
Calibration Function	Yes
Measurement Function	Yes
Unit Dimension	8.0cm (H) x 3.2cm (D)
Unit Weight	117g