

AM73915MZTL Dino-Lite Edge



Dino-Lite Edge 3.0 AM73915MZTL, perangkat USB 3.0 all-in-one dengan optik LWD (Long Working Distance), memungkinkan untuk memenuhi sebagian besar permintaan dan harapan para profesional. Dengan keuntungan dari transmisi kecepatan tinggi pada USB 3.0, AM73915MZTL menghadirkan gambar yang tidak terkompresi dengan 15 FPS pada 5MP atau 45FPS pada resolusi 1.3MP. Ditambahkan dengan fungsi EDOF / EDR, AMR dan FLC canggih yang dapat digunakan dalam aplikasi yang lebih menuntut.

Pertanyaan **contact support** (Klik)

Warranty : 1 Tahun, Made in : Taiwan

Optional :

Basic Stand

Profesional Stand





SPECIFICATIONS


Special feature	Auto Magnification Reading,Ext Depth of Field,FLC (Switch LED)
ESD Safe	Yes
Measurement functionality	Yes
Housing	Aluminum
Cap Changeable	Yes(AD / EDGE series)
Polarizer	Yes
Number of LEDs	8
Light / LED type	White
Working distance	Long
Magnification	10~140x
Resolution	5 MP (2592 x 1944)
Connection	USB 3.0

ADDITIONAL INFO

Specifications

Model	AM73915MZTL Dino-Lite Edge 3.0
Interface	USB 3.0
Product Resolution	5M pixels (2560x1920)
Magnification	10x~140x
Frame Rate	15fps at 5MP, YUY2 (USB 3.0) 45fps at 1.3MP, YUY2 (USB 3.0) 15fps at 1.3MP, YUY2 (USB 2.0)
Lighting	8 white LEDs
Polarizer	Yes

Microtouch	Yes
Operating System Supported	Windows 10, 8, 7 Mac OS 10.9 or later (EDOF/EDR are NOT controllable on Mac OS)
System Requirements	· Intel® Core™ i5 3.2GHz · 20GB available HD space · 8GB available RAM or more · 2GB video card or more · CD-ROM drive for installation · USB 3.0 with intel chipset
Calibration Function	Yes
Measurement Function	Yes
Unit Dimension	11.9cm (H) x 3.3cm (D)
Unit Weight	110g



Informasi tentang jarak kerja dan bidang pandang

M	WD	FOV (x)	FOV (y)	DOF
10	234.5	39.0	29.2	12.0
20	112.5	19.5	14.6	4.5
30	72.5	13.0	9.7	3.1
40	54.5	9.8	7.3	1.8
50	43.5	7.8	5.8	1.2
60	37.5	6.5	4.8	0.9
70	33.7	5.6	4.2	0.6
80	31.3	4.9	3.6	0.5
90	29.9	4.3	3.2	0.37
100	29.2	3.9	2.9	0.28
110	29.1	3.5	2.6	0.22
120	29.3	3.3	2.4	0.17
130	29.8	3.0	2.2	0.12
140	30.5	2.8	2.1	0.09

M = perbesaran tingkat WD = kerja jarak (tanpa tutup depan) FOV = bidang pandang DOF = kedalaman bidang Unit = mm

