



Dino-Lite System Requirements

Required/minimum

- OS: Windows ~~XP, Vista~~, 7, 8, 10 / Mac OS 10.5+*
- CPU: 2GHz Intel Core Series (or equivalent)
- HD: 5GB free
- RAM: 4GB
- Graphics: CPU-embedded "Intel HD" (or equivalent) or discrete 256M video card
- CD-ROM drive or internet connection for software installation
- USB 2.0

Recommended (Disarankan)

- OS: Windows 8, 10, (Win 11 hubungi kami) / Mac OS 10.8+*
- CPU: 2.5 GHz Intel Core i3/i5/i7 Series (or equivalent)
- HD: 8GB free
- RAM: 4GB / 8GB Lebih baik
- Graphics: CPU-embedded "Intel HD" (or equivalent) or discrete 256MB video card / shared memory 2GB
- CD-ROM drive or internet connection for software installation
- **USB 2.0 / (USB 3.0 Wajib untuk series 73xxx)**

*Untuk penggunaan Mac OS harap konsultasi tipe yang sesuai.

Lihat tabel selanjutnya atau [Hubungi kami](tel:08551700051) 08551700051

Update Jan 2022



7025	X	X	✓	✓	✓	✓	✓	✓	✓
7115	✓	✓	✓	✓	✓	✓	✓	✓	✓
7515	X	X	✓	✓	✓	✓	✓	✓	✓
7915	X	X	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡
73115	X	X	✓	✓	✓	✓	✓	✓	✓
73915	X	X	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡

*Partial compatibility/not tested

✓† Partial compatibility, EDOF/EDR are not compatible

✓‡ Partial compatibility, FLC/AMR are compatible. EDOF/EDR are not compatible