



Dino-Lite System Requirements

Required/minimum

- OS: Windows 10 / Mac OS 10.5+*
- CPU: Intel Core i3 series (or equivalent)
- HD: 5GB free
- RAM: 4GB
- Graphics: Intel discrete 2GB video card
- CD-ROM drive or internet connection for software installation
- USB 2.0

**Untuk penggunaan Mac OS harap konsultasi tipe yang sesuai.*

Recommended (Disarankan)

- OS: Windows 10/11 / Mac OS 10.8+*
- CPU: 2.9 GHz dual-core Intel Core i5 (equivalent or greater)
- HD: 10GB free
- RAM: 8GB or more
- Graphics: Intel discrete 4GB video card
- CD-ROM drive or internet connection for software installation
- USB 3.0 & Tipe C

Note : Dino-Lite Series 7 (5 Mega Pixel) dan 8 Mega Pixel harus mengikuti Recommended spec.

Lihat tabel selanjutnya atau [Hubungi kami](#).

Update 2023

support@microscope.id

[0855 1700 051](tel:08551700051)



	10.7	10.8	10.9	10.10	10.11	10.12	10.13	10.14	10.15
2111	✓	✓	✓	✓	✓	✓	✓	✓	✓
3111	✓	✓	✓	✓	✓	✓	✓	✓	✓
3113	✓	✓	✓	✓	✓	✓	✓	✓	✓
3713	X	X	X	X	X	X	X	X	X
3715	X	X	X	X	X	X	X	X	X
4023	✓	✓	✓	✓	✓	✓	✓	✓	✓
4111	✓	✓	✓	✓	✓	✓	✓	✓	✓
4113	✓	✓	✓	✓	✓	✓	✓	✓	✓
4115	✓	✓	✓	✓	✓	✓	✓	✓	✓
4515	X	X	✓	✓	✓	✓	✓	✓	✓
4815	X	X	✓†	✓†	✓†	✓†	✓†	✓†	✓†
4915	X	X	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡

support@microscope.id

[0855 1700 051](tel:08551700051)



7013	✓	✓	✓	✓	✓	✓	✓	✓	✓
7023	✓	✓	✓	✓	✓	✓	✓	✓	✓
7025	X	X	✓	✓	✓	✓	✓	✓	✓
7115	✓	✓	✓	✓	✓	✓	✓	✓	✓
7515	X	X	✓	✓	✓	✓	✓	✓	✓
7915	X	X	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡
73115	X	X	✓	✓	✓	✓	✓	✓	✓
73915	X	X	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡	✓‡

*Partial compatibility/not tested

✓† Partial compatibility, EDOF/EDR are not compatible

✓‡ Partial compatibility, FLC/AMR are compatible. EDOF/EDR are not compatible

support@microscope.id

[0855 1700 051](tel:08551700051)