

HXD-V OND

## Backlight for Microscope

### Microscope Backlight Light Source

- Adjustable brightness.
- The outside diameter of ring light: Ø95mm.
- The inner diameter of ring light: Ø72mm
- Color temperature 6000K



The extremely LED Backlight is easy to use and fits most bases on the market. Backlighting is one of the most popular forms of vision illumination for microscopes, providing maximum contrast between object and scene and showing a high contrast image with consistent edge definition

- Item Type: Microscope Bottom Light Source HXD-V OND

- Input Voltage: AC 110V-240V

- Output Voltage: 24V (Max)

- Output Power: 9W

- Inner Diameter: 70 mm

- Outer Diameter: 95 mm

- Height: 18 mm

- Light model : LED

- Brightness Adjustable

- Working Environment: temperature: 0-40 degree, humidity: 20-85% RH

- Suitable for: **stereo microscope, monocular microscope and other automation equipment**

Made in China



Complex Pertokoan Glodok Plaza Blok F No. 115, Lt4 Jakarta 11180 PO. Box 4978 JKT / 11049 Indonesia  
Telp : 0855-1700-051

Mail : [Support@Microscope.id](mailto:Support@Microscope.id) Website : [www.microscope.id](http://www.microscope.id)