


RING LIGHT

for microscope

microscope.id



Examples :



FEATURES :

- ✓ 144 White LED bulbs
- ✓ 4 Zones with individual ON/OFF
- ✓ 10 Adjustable brightness control
- ✓ Strong and focused illumination with shadow free

SPECIFICATIONS :

- ✓ 144 White LED bulbs
- ✓ 4 Zones with individual ON/OFF
- ✓ 10 adjustable mode for brightness control
- ✓ Outer diameter: 98 mm (3.86 inch)
- ✓ Inner diameter: 61 mm (2.40 inch)
- ✓ Power: 4.5W
- ✓ Input Voltage: 90V - 265V AC
- ✓ Output Voltage: 12V DC
- ✓ Brightness Control: 0 ~ 100%
- ✓ Working Distance: 40 ~ 250 mm
- ✓ Brightness Level: 200,000 Lux. from 100 mm high
- ✓ Color Temperature: 6400K
- ✓ Supplied with an additional 48 mm (inner diameter) ring adapter
- ✓ Dimension of main panel box: approx. 60 (L) x 110 (W) x 32 (D) mm (2.36 x 4.33 x 1.26 inch)